



URGO AG



URGO CLEAN



URGO START



URGO TOL



Dressing Cutting Guide

Contents

Key	p3
Benchmarks in wound healing	p4
Technology Lipido-Colloid	p5
Optimize the wound bed to initiate healing	p6
Tailor the treatment to each patient	p7
Hands & fingers	
• Hand	p8
• Whole hand	p9
• Fingers	p10
• Between fingers	p11-12
• Palm	p13
• Thumb	p14
Knees & feet	
• Diabetic foot	p15
• Knee	p16
• Toes	p17
• Heel	p18 - 19
• Ankle	p20
Face & head	
• Upper lip	p21
• Nose	p22
• Ear	p23
• Chin	p24
Elbows & underarm	
• Elbow	p25 - 26
• Underarm	p27 - 28
Tubing/wound drainage sites	p29
Products	p30 - 31
References	p32

Key



Hands and Fingers



Non to low exuding wounds



Knees



Low to moderately exuding wounds



Feet



Moderately to heavily exuding wounds



Heel



Face and Head



Elbow



Underarm



Tubing

Benchmarks in wound healing

Wound healing is a complex phenomenon and each phase poses particular problems.

In two stages and four steps, Urgo Medical offer powerful and innovative solutions adapted to each type of wound and patient :



OPTIMIZE THE WOUND BED TO INITIATE HEALING



In case of risk of infection, obtain a safe large-spectrum bactericidal effect



UrgoTul Ag/Silver - UrgoCell Ag/Silver



Achieve efficient and atraumatic desloughing



UrgoClean



TAILOR THE TREATMENT TO EACH PATIENT



Accelerate wound healing



UrgoStart



Optimise wound healing to achieve the best possible aesthetic and functional results




UrgoTul

Technology Lipido-Colloid



TLC healing matrix

Composition of the TLC



CMC particles

Lipophilic particles

TLC is the core technology of all the Urgo Medical products. It is indeed compatible with different brackets (interface, fibers ...) and can be associated with different compounds (NOSF, Ag ...)

Benefits of the TLC healing matrix

- **Allows a solid and harmonious reconstruction of the dermis and epidermis**
 - Stimulates fibroblasts and extra-cellular matrix components proliferation
 - Moist wound healing
- **Allows pain-free and atraumatic care**
 - Keeps its greasy enduction over time, does not dry out
 - Does not adhere to the wound
 - High tolerance



Combination of TLC and NOSF: a unique compound designed to neutralise excess MMPs: **to rebalance wound equilibrium and promote faster healing**



Combination of TLC and Silver: antibacterial agent: **to manage infection.**

Optimize the wound bed to initiate healing




In case of risk of infection, obtain a safe large-spectrum bactericidal effect



Thin absorbent dressing with TLC-Ag



Thin absorbent dressing with TLC-Ag and adhesive border



Foam dressing with TLC-Ag



Foam dressing with TLC-Ag and adhesive border



Achieve efficient and atraumatic desloughing



Hydro-desloughing soft-adherent absorbent dressing with TLC

Tailor the treatment to each patient



Accelerate wound healing



Foam dressing with soft-adherent TLC-NOSF



Flexible contact layer with TLC-NOSF



Foam dressing with soft-adherent TLC-NOSF and adhesive border



Optimise wound healing to achieve the best possible aesthetic and functional results



Thin foam dressing with soft-adherent TLC



Thin absorbent dressing with TLC and adhesive border



Foam dressing with soft-adherent TLC



Foam dressing with soft-adherent TLC and adhesive border



Hand




PREPARE



UrgoTul Ag Life



UrgoClean



HEAL



UrgoStart Contact



UrgoTul Life

- Fold the dressing in half lengthwise and cut out the pattern marked by the dotted lines on the diagram.
- Apply the dressing to the wound.
- Cover with a secondary dressing for UrgoStart Contact.
- Secure the dressing with tape or a retention bandage.



Whole hand




PREPARE



UrgoTul Ag Life



UrgoClean



HEAL



UrgoStart Contact



UrgoTul Life

- Fold the dressing in half lengthwise and cut out the pattern marked by the dotted lines on the diagram.
- Apply the dressing to the wound.
- Cover with a secondary dressing for UrgoStart Contact.
- Secure the dressing with tape or a retention bandage.



Fingers




PREPARE



UrgoTul Ag Life



UrgoClean



HEAL



UrgoStart Contact



UrgoTul Life

- Cut out the pattern marked by the dotted lines on the diagram.
- Apply the dressing to the wound.
- Cover with a secondary dressing for UrgoStart Contact.
- Secure the dressing with tape, film dressing or a retention bandage.



Between fingers




PREPARE



UrgoTul Ag Life



UrgoClean



HEAL



UrgoStart Contact



UrgoTul Life


- Fold the dressing in half lengthwise and cut out the pattern marked by the dotted lines on the diagram.
- Adapt the pattern to the space between the fingers.
- Apply the dressing over the hand.
- Cover with a secondary dressing for UrgoStart Contact.
- Secure with a retention bandage as appropriate.



Between fingers



PREPARE



HEAL



UrgoStart (Heel)



UrgoTul Absorb (Heel)

- Cut out the pattern marked by the dotted lines on the diagram.
- Adapt the pattern to the space between the fingers and then carefully peel off the protective backing without touching the dressing.
- Apply the dressing over the hand, starting with the fingertips.
- Secure the dressing with a retention bandage as appropriate.



Palm

PREPARE



UrgoTul Ag Life Border 

UrgoCell Ag Border 



HEAL



UrgoStart Border 



UrgoTul Life Border 


UrgoTul Absorb Border 

- In order to ensure a good bond of the dressing to the base of the fingers cut the edges first, then apply.



Thumb

PREPARE



HEAL



- Cut the dressing into a T shape as shown and position the thumb perpendicularly to the dressing.
- Pull back each side of the edges, in the number order shown on the diagram, to cover the whole thumb.




Foot (diabetic foot ulcers)



HEAL



- 1 : Under the toe (pulp)
- 2 : Under the foot
- 3 : Between toes
- 4 : Dorsum of the foot
- 5 : Internal side of the foot
- 6 : External side of the foot
- 7 : Heel



- Prepare the wound
 1. Gently remove the slough if present.
 2. Clean the wound with normal saline.
 3. If an antiseptic is used first, rinse the wound thoroughly with saline solution before applying UrgoStart Contact.
- Apply the dressing
 1. Flexible and very conformable, UrgoStart Contact is especially suitable for lining cavity or awkwardly placed wounds and for wicking deep wounds. UrgoStart Contact can be cut using sterile equipment to fit the dressing to the wound.
 2. Cover UrgoStart Contact with a secondary dressing suitable to absorb the level of wound exudate and hold the secondary dressing in place with tape or a retention bandage.
- Offload pressure areas



Knee



PREPARE



UrgoTul Ag/Silver



UrgoClean

HEAL




UrgoStart



UrgoTul Absorb

- Cut out the pattern marked by the dotted lines on the diagram.
- Apply the dressing over the knee.
- Secure the dressing with a retention bandage as appropriate.



Toes



PREPARE



UrgoCell Ag/Silver



UrgoClean

HEAL




UrgoStart



UrgoTul Absorb


- Cut out the pattern marked by the dotted lines on the diagram.
- Apply the dressing over the toe.
- Secure the dressing with a retention bandage as appropriate.



Heel

PREPARE

 UrgoTul Ag Lite Border 
 UrgoCell Ag Border 



HEAL

 UrgoStart Border 

 UrgoTul Lite Border 
 UrgoTul Absorb Border 




- Cut out a V shape on both sides, slightly overlapping the pad.
- This will aid the dressing to shape to the heel without folding.

Heel



PREPARE



HEAL

 UrgoStart (Heel)

 UrgoTul Absorb (Heel)




- Apply the dressing over the heel.
- Secure the dressing with a retention bandage as appropriate.



Ankle

PREPARE



HEAL




- Cut out a V shape on both sides, slightly overlapping the pad.
- This will aid the dressing to shape to the heel without folding.



Upper lip



PREPARE



HEAL



- Fold the dressing in half lengthwise and cut out the pattern marked by the dotted lines on the diagram.
- Adapt the cut-out to the contour and size of upper lip.
- Secure the dressing with tape.



Nose




PREPARE



UrgoTul Ag Life



UrgoClean



HEAL



UrgoStart Contact



UrgoTul Life

- Fold the dressing in half lengthwise and cut out the pattern marked by the dotted lines on the diagram.
- Adapt the outline to the shape of the nose.
- Cover with a secondary dressing for UrgoStart Contact.
- Secure the dressing with tape.



Ear




PREPARE



UrgoCell Ag Life



UrgoClean



HEAL



UrgoStart Contact



UrgoTul Life

- Fold the dressing in half lengthwise and cut a hole into the fold by following the dotted line on the diagram.
- Position the ear inside the hole and bend the fold upwards.
- Cover with a secondary dressing for UrgoStart Contact.
- Secure the dressing with a net dressing.



Chin




PREPARE



UrgoTul Ag Lite



UrgoClean



HEAL



UrgoStart Contact



UrgoTul Lite

- Fold the dressing in half lengthwise. The cutting shape depends on the wound surface area.
- Cover with a secondary dressing for UrgoStart Contact.
- Secure the dressing with tape.



Elbow


PREPARE



UrgoCell Ag Lite Border



UrgoCell Ag Border



HEAL



UrgoStart Border



UrgoTul Lite Border



UrgoTul Absorb Border




- Cut out the V shape in the adhesive border to shape the elbow.
- Apply dressing when the arm is sem-bent to avoid stretching when moving.



Elbow



PREPARE



HEAL



UrgoStart (Heel)



UrgoTul Absorb (Heel)

- It is not necessary to cut the dressing to size.
- Carefully peel off the protective backing without touching the dressing.
- Apply the dressing so that the tip points towards the hand.
- Secure the dressing with a retention bandage.



Underarm

PREPARE



UrgoCell Ag Lite Border



UrgoCell Ag Border



HEAL





UrgoStart Border



UrgoTul Lite Border



UrgoTul Absorb Border



- Cut out the corners of the adhesive border.
- Apply as shown to ensure ease of movement.



Underarm




PREPARE



UrgoCell Ag/Silver



UrgoClean



HEAL



UrgoStart



UrgoTul Absorb

- Cut out the pattern marked by the dotted lines on the diagram.
- Carefully peel off the protective backing without touching the dressing.
- Apply the dressing to the wound and secure it in place with tape or a retention bandage as appropriate.



Tubing/wound drainage sites





PREPARE



UrgoCell Ag/Silver



UrgoClean



HEAL



UrgoStart



UrgoTul Absorb

- Cut out the pattern marked by the dotted lines on the diagram.
- Carefully peel off the protective backing without touching the dressing.
- Apply the dressing and then secure it with a film dressing or tape.

Products



	Sizes	Absorption
UrgoTul Ag Life	10x12 cm 15x20 cm	1
UrgoTul Ag Life Border	6.5x10 cm 8x8 cm 10x12 cm	1
UrgoCell Ag/Silver	6x6 cm 10x10 cm 10x12 cm 15x15 cm 15x20 cm	2
UrgoCell Ag Border	13x13 cm	2



	Sizes	Absorption
UrgoClean	6x6 cm 10x10 cm 10x12 cm 15x15 cm 15x20 cm	1



	Sizes	Absorption
UrgoStart	6x6 cm 10x10 cm 10x12 cm 13x12 cm 15x15 cm 15x20 cm Heel 12x19 cm	2
UrgoStart Contact	5x7 cm 10x10 cm 10x12 cm 13x12 cm 15x20 cm	1



	Sizes	Absorption
UrgoStart Border*	8x8 cm 10x10 cm 12x12 cm 13x13 cm 15x15 cm 15x20 cm Sacrum 20x20 cm	3



	Sizes	Absorption
UrgoTul Life	6.5x10 cm 13x12 cm 10x20 cm 15x20 cm	1
UrgoTul Life Border	6.5x10 cm 8x8 cm 10x12 cm 10x20 cm 15x20 cm	1
UrgoTul Absorb	6x6 cm 10x12 cm 13x12 cm 15x15 cm 15x20 cm Heel 12x19 cm	2
UrgoTul Absorb Border*	8x8 cm 10x10 cm 12x12 cm 13x13 cm 15x15 cm 15x20 cm Sacrum 20x20 cm	3

* The adhesive edges can be cut. However, as it includes a super-absorbent layer, the central pad should not be cut.

SKIN Wounds

Innovation


Listening Evidences

SOLUTIONS

Expertise

Commitment

Technology



Our teams are impassioned by wound healing research
because they are inspired by life.



TO HEAL

To live TO REPAIR

TO LOVE

To erase To touch

TO MOVE

TO SHARE

To cure

Today, being comfortable in one's own skin is vital to well-being. The physical and psychological effects of wounds should not interfere in the patient's ability to move, act and feel. Knowing this, Urgo Medical has developed a unique approach to wound healing research:

- **Leading-edge knowledge of skin and its mechanisms:** every day, our researchers bring progress to tissue repair;
- **A vision incorporating all dimensions of healing:** therapeutic, practical, functional, aesthetic, psychological, etc.

That is why we are so close to patients and care providers. They inspire us every day.

 **URGO**
MEDICAL
HEALING PEOPLE®